

**NAME**

pyFAI-calib – Geometry calibration tool of PyFAI

**SYNOPSIS**

**pyFAI-calib** [*options*] [*images*]

**DESCRIPTION**

pyFAI-calib

A tool for determining the geometry of a detector using a reference sample.

**OPTIONS**

- h, --help**  
show this help message and exit
- V, --version**  
print version of the program and quit
- o FILE, --out=FILE**  
Filename where processed image is saved
- v, --verbose**  
switch to debug/verbose mode
- g GAUSSIAN, --gaussian=GAUSSIAN**  
Size of the gaussian kernel. Size of the gap (in pixels) between two consecutive rings, by default 100 Increase the value if the arc is not complete; decrease the value if arcs are mixed together.
- c, --square**  
Use square kernel shape for neighbour search instead of diamond shape
- b BACKGROUND, --background=BACKGROUND**  
Automatic background subtraction if no value are provided
- d DARK, --dark=DARK**  
list of dark images to average and subtract
- f FLAT, --flat=FLAT**  
list of flat images to average and divide
- r, --reconstruct**  
Reconstruct image where data are masked or <0 (for Pilatus detectors or detectors with modules).
- s SPLINE, --spline=SPLINE**  
spline file describing the detector distortion
- p PIXEL, --pixel=PIXEL**  
size of the pixel in micron
- D DETECTOR\_NAME, --detector=DETECTOR\_NAME**  
Detector name (instead of pixel size+spline)
- m MASK, --mask=MASK**  
file containing the mask (for image reconstruction)
- n NPT, --npt=NPT**  
file with datapoints saved