

**NAME**

pyFAI-recalib – Geometry refinement / recalibration tool of PyFAI

**SYNOPSIS**

**pyFAI-recalib** [*options*] [*images*]

**DESCRIPTION**

pyFAI-recalib

An automatic tool for refining the geometry of a detector using a reference sample and a good guess.

**OPTIONS**

- h, --help**  
show this help message and exit
- V, --version**  
print version of the program and quit
- o FILE, --out=FILE**  
Filename where processed image is saved
- v, --verbose**  
switch to debug mode
- s FILE, --spacing=FILE**  
file containing d-spacing of the reference sample
- d FILE, --dark=FILE**  
list of dark images to average and subtract
- f FILE, --flat=FILE**  
list of flat images to average and divide
- m FILE, --mask=FILE**  
file containing the mask
- p FILE, --poni=FILE**  
file containing the diffraction parameter (poni-file)
- n FILE, --npt=FILE**  
file with datapoints saved
- e ENERGY, --energy=ENERGY**  
energy of the X-Ray beam in keV (hc=12.398keV.A)
- w WAVELENGTH, --wavelength=WAVELENGTH**  
wavelength of the X-Ray beam in Angstrom
- l DISTANCE, --distance=DISTANCE**  
sample-detector distance in millimeter
- poni1=PONI1**  
poni1 coordinate in meter
- poni2=PONI2**  
poni2 coordinate in meter
- rot1=ROT1**  
rot1 in radians
- rot2=ROT2**  
rot2 in radians
- rot3=ROT3**  
rot3 in radians